



DTX Studio™ Clinic

Version 4.10

Silent Installation

Table of Contents

Scope	3
Silent Installation Instructions (Windows)	4
InstallShield Instructions	4
Windows Device Drivers	4
Skip Device Driver Installation	5
Devices Supported by the Windows Device Drivers	5
Silent Installation Instructions for Mac?	5
Silent Configuration	6
Prerequisites	6
Clean Installation	6
Upgrading DTX Studio Clinic	6
Command Line Interface	7
Usage	7
Installation Key	8
Get the Installation Key	8
Data Directory	8
Force Flag	8
EULA	8
Settings File	9
DTX Studio Core Connection	9
Examples	9
Full Configuration	9
Updating an Existing Installation to a New License	10
Setting the User Settings Only	10
Set up DTX Studio Core Connection with the Specified URL	10
Set up DTX Studio Core Connection with the Core URL as Configured on DTX Studio Go	10
Clear the DTX Studio Core Connection Configuration	10
Only accept the EULA (e.g. after the upgrade of an existing installation)	10
Output	11
SCCM Application Creation Guide (Windows)	13
Prerequisites	13
Single Application Package with PowerShell Script	13
Technical Support	19
Annex I: End User License Agreement (EULA)	20


Scope

This procedure was created to allow enterprise customers to install DTX Studio Clinic on a number of machines from a central point of control. This also facilitates batch installations and upgrades via a remote administrative tool, like Microsoft System Center Configuration Manager (SCCM).

This document sets apart the administrative aspect (operating system specific installation procedures) and functional aspect post-installation (DTX Studio Clinic configuration). It attempts to create a comprehensive overview of the general installation procedure for a large practice with centralized administration workflow.

In short, the silent installation procedure is a two-step approach to be done by a system administrator:

- Install using InstallShield command line (**silent installation**). This is currently only supported for Windows networks.
- Run the first time via command line and with administrative privileges to do the **silent configuration** procedure. This step is OS independent.



Silent Installation Instructions (Windows)

The silent installation done by a system administrator on a Windows network typically consists out of the following steps:

- Configure registry keys to prevent device driver packages being installed.
- Install DTX Studio Clinic using the InstallShield command line (silent installation).
- Run the software for the first time via command line and with administrative privileges to perform the silent configuration.
- Install specific device driver(s) needed on workstations that need to support USB devices.

InstallShield Instructions

To launch the installer in silent mode, start it with the following command line parameters:

```
DTXStudioClinic(x64)_x.x.x.x.exe /s /v"/qn"
```

Normally, this command will not trigger a machine restart. However, if you want to guarantee that the computer does not reboot, this behaviour can be explicitly enforced:

```
DTXStudioClinic(x64)_x.x.x.x.exe /s /v"/qn"REBOOT=ReallySuppress
```

Windows Device Drivers

DTX Studio Clinic attempts to install Windows device drivers for supported USB intraoral cameras during installation. Driver installation can be disabled during rollout if device drivers are not required on all workstations. For computers that do require device drivers, either the installer must be run with driver installation enabled, or individual driver packages can be installed manually as a post-installation step.

DEXIS Software Manager, which can optionally be installed alongside DTX Studio Clinic, is the preferred method for deploying USB intraoral sensor drivers. However, it is not included as part of the silent installation and must be installed separately. The installer does not automatically install USB intraoral sensor drivers during the initial setup. Instead, users are notified when a supported sensor is connected, and the drivers can be installed either in advance or afterwards on the machines that require them. Individual driver packages are available on request, and DEXIS intraoral sensor drivers can also be downloaded from the DEXIS Download Center. All drivers are digitally signed by either Microsoft or the device manufacturer.

The DTX Studio Camera Drivers prerequisite package installs the required certificates, preventing users from being prompted to approve driver installations during deployment.

If driver deployment issues occur, deployment scripts and driver certificates can be provided separately. Alternatively, problematic driver certificates can be imported into the Trusted Publishers certificate store before installing the drivers, which prevents additional installation prompts.

For environments where drivers are deployed manually, individual driver packages are available on request.

- **DTX Studio Camera Drivers(x64)_x.x.x.x.exe** (for the intraoral camera drivers)
- **DEXIS_Intraoral_Sensor_Driver_Setup_Launcher_x.x.x.x.exe** (for the intraoral sensor drivers via the DEXIS Software Manager).

Skip Device Driver Installation

If you do not need the device drivers on all client machines, the device driver installation can be skipped. This can be achieved by setting the following registry value (DWORD) to **1** before installing DTX Studio Clinic:

HKEY_LOCAL_MACHINE\SOFTWARE\Medicim\DTX Studio Clinic\SkipCameraDrivers

In a script, use the following command:

```
REG ADD "HKEY_LOCAL_MACHINE\SOFTWARE\Medicim\DTX Studio Clinic"/v  
SkipCameraDrivers /t REG_DWORD /d 1 /f
```

Devices Supported by the Windows Device Drivers

– Intraoral sensors:

DEXIS™ Ti2, DEXIS™ Ti2- S, DEXIS™ Titanium, DEXIS™ IXS,
DEXIS™ IXS - S, DEXIS™ Platinum, GXS-700™, Gendex™ GXS-700™

– Intraoral cameras:

DEXIS™ DexCAM™ 4 HD, DEXIS™ DexCAM™ 3, DEXIS™ DexCAM™ 4, Gendex™ GXC-300™,
KaVo DIAGNOcam Vision Full HD, KaVo DIAGNOcam™, CariVu™, KaVo ERGOcam™ One.

Silent Installation Instructions for Mac?

Currently there is no solution available for a silent installation of the Mac .pkg installer. If this is a requirement for your company, please contact DTX Studio Technical Support (support@dtxstudio.com).

Silent Configuration

DTX Setup is the application to configure DTX Studio Clinic. The DTX Setup executable can be run in silent configuration mode.

Important

Silent installation is not supported for offline activated licenses.

Prerequisites

An internet connection is needed to successfully validate an installation of the DTX Studio Clinic software.

Make sure to always close all instances of DTX Studio Clinic before attempting this silent configuration.

Running the executable with the silent configuration parameters while DTX Studio Clinic is still running, will not work. Typically, though, the configuration does not need to be executed again when updating an existing installation, so in most cases this is no issue.

Clean Installation

After a clean installation, invoking the DTX Studio Clinic executable with specific command line parameters will launch the application in silent configuration mode, running the setup wizard without making use of the GUI. This is only needed for a fresh install (See '[Upgrading DTX Studio Clinic](#)' on page 6).

Calling the executable with these parameters must be done with administrative privileges.

Upgrading DTX Studio Clinic

In an upgrade scenario, the existing license can be used with the new version, if still valid. However, the new EULA needs to be formally accepted with **-e EULA=ACCEPT** before the first run after upgrading. The only exceptions might be wanting to deploy a new shared settings file (see '[Settings File](#)' on page 9) or changing to another license (see '[Force Flag](#)' on page 8).

Command Line Interface

There is a set of command line parameters that allow users/administrators to configure their installation of DTX Studio Clinic:

Parameter	Description	Information
-k	Installation key	Installation key from DTX Studio License Server (See ' Installation Key ' on page 8).
-d	Data directory	Data location for patient data and program data (See ' Data Directory ' on page 8).
-f	Force installation	Overwrite existing valid license (See ' Force Flag ' on page 8).
-e	Accept EULA	Parameter to accept the End User License Agreement (EULA) (See ' EULA ' on page 8).
-s	Settings file	Shared settings file (See ' Settings File ' on page 9).
-t	Core URL	URL for DTX Studio Core (See ' DTX Studio Core Connection ' on page 9).
-remove	Remove device configuration	Remove device configuration from settings. Used for TWAIN devices, file watcher folders and locally disabled devices.

Usage

When configuring the installation, **-k** and **-e** are mandatory and **-d**, **-s**, **-t**, are optional.

The **-k** argument always needs to be accompanied with the **-e** flag, as the EULA always needs to be formally accepted when taking a license. Consequently, **-d** and **-e** are only possible with the **-k** parameter (installation).

When switching to another license, hence with a different install key - use the **-f** parameter to force a switch of licensing.

The user settings can be passed with **-s** (path to .set file) on initial install. Or as a stand-alone command with a single **-s** as parameter (provided that you already have a valid installation / license).

When upgrading to a newer version of the software, the silent configuration does not need to run again, provided that the previous parameters still apply (e.g. same license, same data directory, same credentials, ...).

Installation Key

The installation key is an alpha-numeric string (e.g. AAB3D-3ABAE-5D4AF-4929A12B14-8771C05), checked by the license server, which allows you to register your installation online. If there is already a valid software license present, the provided installation key will not be used unless the **-f** flag is present to force it (See '[Force Flag](#)' on page 8).

Get the Installation Key

To get the key, go to DTX Studio Go (go.dtxstudio.com). The portal is best viewed with Chrome or Edge browsers. On the [Software](#) page, click the DTX Studio Clinic license. On the DTX Studio Clinic license page, copy the installation key.

Data Directory

The installation data directory is the location of a single central directory where DTX Studio Clinic saves all of its application data. This should be the full path to a location on a local disk (not a network folder!), accessible to all users who might log on to the machine. All data in this folder will be stored in an encrypted way. The folder will be created if it does not exist yet and user permissions will be set accordingly. If nothing provided, the default path will be applied:

C:\ProgramData\DTX Studio\Clinic on Windows

/Users/Shared/DTX Studio/Clinic/ on MacOS.

Force Flag

This flag indicates that the application should override the previous license, even if it is still valid. This is mainly used for a license change (switch to different license with other installation key).

EULA

The EULA must be accepted, both in a clean installation and an upgrade. You must read the EULA (See '[Annex I: End User License Agreement \(EULA\)](#)' on page 20) and include **-e EULA=ACCEPT**.

By including this parameter you confirm the acceptance of the licence by you and your organization. In absence of the inclusion of the EULA parameter the installation or upgrade fails.

Note

The IT administration managing the installation of the software remotely for their users must notify them of the contents of the End User License Agreement.

Settings File

The settings file is the path to a DTX Studio Clinic shared settings file (.set), which can be exported by another DTX Studio Clinic installation. This file holds the user settings (preferences) to be applied to the new installation, allowing for an easy roll-out of user preferences to all the DTX Studio Clinic installations within a practice. If not available, the default settings will be used.

DTX Studio Core Connection

When passing the `-t` parameter, the DTX Studio Core connection is automatically set up. The connection may have to be approved for each workstation on the DTX Studio Core UI page.

If `-t` is not specified, the installation will fall back to the URL that is associated with your DTX Studio Core license on the license server, if it is configured.

There are different possible values for the parameter value:

- A valid URL (e.g. "`http://coreserver:26850`")
- `""`: use the DTX Studio Core URL as configured on DTX Studio Go.
- `"clear"`: clear the DTX Studio Core configuration.

If the passed-on URL or the URL on DTX Studio Go is invalid, error 19 will be returned.

Valid URLs to the DTX Studio Core server are with `http` or `https` scheme, and with or without endpoint specified.

`https://coreserver:26850` and `https://coreserver:26850/dw-endpoint/api/`, for example, are both valid URLs.

Examples

Full Configuration

Run the following command, for example, with administrative privileges:

```
DTXStudio.exe -k "AAB3D-3ABAE-5D4AF-4929A-12B14-8771C05" -f -e EULA=ACCEPT
-d "C:\ProgramData\DTX Studio\Clinic" -s "C:\Users\John.Doe\Documents\
DTXStudioSettings.set" -t "http://coreserver:26851"
```

This will result in an installation and configuration of DTX Studio Clinic.

The license key used for registration is "AAB3D-3ABAE5D4AF-4929A-12B14-8771C05".

The EULA is accepted.

Its data directory will be set to "C:\ProgramData\DTX Studio\Clinic".

Settings will be loaded from "C:\Users\John.Doe\Documents\DTXStudioSettings.set".

It will connect to "http://coreserver:26851", while overriding the already present installation on that machine, if one exists.

Updating an Existing Installation to a New License

```
DTXStudio.exe -k "340B7-634F6-CA46C-29021-5F0EE-2AAF9AC" -e EULA=ACCEPT -f
```

Setting the User Settings Only

```
DTXStudio.exe -s "C:\Users\John.Doe\Documents\DTXStudioSettings.set"
```

Set up DTX Studio Core Connection with the Specified URL

```
DTXStudio.exe -t "http://coreserver:26851"
```

Set up DTX Studio Core Connection with the Core URL as Configured on DTX Studio Go

```
DTXStudio.exe -t ""
```

Clear the DTX Studio Core Connection Configuration

```
DTXStudio.exe -t "clear"
```

Only accept the EULA (e.g. after the upgrade of an existing installation)

```
DTXSetup.exe -e EULA=ACCEPT
```

Output

The configuration procedure provides a return code to the caller, providing more information on the outcome of the installation.

0	Success
1	Invalid program arguments
2	Failed to create data directory
3	Invalid file system permissions
4	Insufficient administrative privileges
5	Failed to save license
6	Application already running
7	EULA was not formally accepted
8	No connection to DTX Studio License Server
9	License does not exist
10	Maximum number of installations exceeded on license
11	Invalid software release on license
12	License has been disabled
13	License period has expired
14	License has no associated users
15	Valid license already present
16	Settings (file) could not be imported
17	Settings file is empty
18	Imported settings could not be saved
19	Invalid DTX Studio Core URL provided (as provided via -t parameter or as configured on DTX Studio Go). Expected is a valid URL according to the URI standard and supported scheme (http or https).
20	No connection to DTX Studio Core could be made (to the url provided via -t or as configured on DTX Studio Go)
22	Unspecified error
23	No connection to DTX Studio Core could be made while connecting with https and there's no valid certificate installed.
24	No connection to DTX Studio Core could be made as the application is already registered but the authentication info is missing. To solve this, delete this workstation on the DTX Studio Core Workstations page and try again.
255	Internal error (error while running DTX Studio Clinic from DTXSetup, possible only when at least one of these parameters is used: -s, -t, --remove).

When a silent installation results in a return code **3**, it means that either access to the registry or plist has failed, or access to the folder specified as the data directory was not granted.

Output

The specific result code **5** indicates that we were unable to save the result from activation and validation of a license to the registry or plist file.

The specific result code **16** will be returned when a general error occurs during importing of a settings file (.set), or when attempting to import whilst not having a valid installation.

SCCM Application Creation Guide (Windows)

This step-by-step guide helps configuring DTX Studio Clinic for SCCM (Microsoft System Center Configuration Manager).

Prerequisites

The DTX Studio Clinic installer has dependencies on Visual Studio redistributables. It will install (or upgrade) them if not yet present on the system, but this may trigger unwanted reboots during the SCCM installation. For this reason we recommend creating SCCM packages for these dependencies and set them as requirements for the DTX Studio Clinic SCCM package.


Currently the dependencies are:

- Microsoft Visual C++ 2012 Update 4 Redistributable (x64)
- Microsoft Visual C++ 2022 Redistributable (x64)
- Microsoft Visual C++ 2022 Redistributable (x86)

You can get the latest versions from the official Microsoft website.

Single Application Package with PowerShell Script

1. Copy the source files to a network location.
2. Click the [Create Application](#) icon.
3. Select [Manually specify the application information](#).



4. Fill in the **General Information** details.

The screenshot shows the 'Create Application Wizard' window, specifically the 'General Information' step. The left sidebar lists the steps: General, General Information (selected), Software Center, Deployment Types, Summary, Progress, and Completion. The main area is titled 'Specify information about this application'. It contains the following fields and controls:

- Name:** DTX Studio Clinic
- Administrator comments:** (empty text box)
- Publisher:** Medicum NV
- Software version:** 1.8.5.2
- Optional reference:** (empty text box)
- Administrative categories:** (dropdown menu with a 'Select...' button)
- Date published:** ☐ 9/30/2020
- Specify the administrative users who are responsible for this application:**
 - Owners:** eb.sys (with a 'Browse...' button)
 - Support contacts:** eb.sys (with a 'Browse...' button)

At the bottom, there are navigation buttons: '< Previous', 'Next >', 'Summary', and 'Cancel'.

5. Fill in the **Software Center** details.

The screenshot shows the 'Create Application Wizard' window, specifically the 'Software Center' step. The left sidebar lists the steps: General, General Information, Software Center (selected), Deployment Types, Summary, Progress, and Completion. The main area is titled 'Specify the Software Center entry'. It contains the following fields and controls:

- Selected language:** English (United States) default (with an 'Add/Remove...' button)
- Localized application name:** DTX Studio Clinic
- User categories:** (empty text box with an 'Edit...' button)
- User documentation:** (empty text box with a 'Browse...' button)
- Link text:** (empty text box)
- Privacy URL:** (empty text box)
- Localized description:** (empty text box)
- Keywords:** (empty text box)
- Icon:** (a small icon with a 'Browse...' button)
- ☐ Display this as a featured app and highlight it in the company portal

At the bottom, there are navigation buttons: '< Previous', 'Next >', 'Summary', and 'Cancel'.

6. Set the **Deployment Type** to **Script Installer**.

The screenshot shows the 'Create Deployment Type Wizard' dialog box with the 'General' tab selected. The left sidebar lists the steps: General, General Information, Content, Detection Method, User Experience, Requirements, Dependencies, Summary, Progress, and Completion. The main area is titled 'Specify settings for this deployment type'. It contains a 'Type:' dropdown menu set to 'Script Installer'. Below it, there are two radio buttons: 'Automatically identify information about this deployment type from installation files' (unselected) and 'Manually specify the deployment type information' (selected). A 'Location:' text box is present with a 'Browse...' button and an example path: '\\Server\Share\File'. At the bottom, there are four buttons: '< Previous', 'Next >', 'Summary', and 'Cancel'.

7. Enter the details of the **Deployment General Information**.

The screenshot shows the 'Create Deployment Type Wizard' dialog box with the 'General Information' tab selected. The left sidebar lists the steps: General, General Information, Content, Detection Method, User Experience, Requirements, Dependencies, Summary, Progress, and Completion. The main area is titled 'Specify general information for this deployment type'. It contains a text box for 'Name' with the value 'Installation'. Below it is a text box for 'Administrator comments'. At the bottom, there is a 'Languages:' section with a list box and a 'Select...' button. At the bottom of the dialog, there are four buttons: '< Previous', 'Next >', 'Summary', and 'Cancel'.


8. Enter the full UNC path to the network location of the source files for the Content Location. If you used the sample PowerShell script, the Installation Program is
- ```
powershell.exe -executionpolicy bypass -nologo -nopprofile -file ".\Install_DTX_Studio.ps1"
```

The Uninstall Program is


```
MsiExec.exe /X{816FA862-BBC7-41EC-A48B-AE4AE6C3226C} /qn /norestart
```

**Note**

The product ID may change depending on the DTX Studio Clinic version being installed. Check the registry for the correct product ID.



9. Configure the detection rule [Setting Type](#) to [Windows Installer](#). Enter product ID from the registry. It will be the same Product ID used by the uninstall string in the previous step.



10. Configure the **User Experience** settings.

Specify user experience settings for the application

Installation behavior: Install for system if resource is device; otherwise install for user

Logon requirement: Whether or not a user is logged on

Installation program visibility: Hidden

☐ Allow users to view and interact with the program installation

Specify the maximum run time and estimated installation time of the deployment program for this application. The estimated installation time displays to the user when the application installs.

Maximum allowed run time (minutes): 120

Estimated installation time (minutes): 10

< Previous   Next >   Summary   Cancel

11. Configure the **Requirements** to **Windows 10 64-bit** or **Windows 11 64-bit**.

Specify installation requirements for this deployment type

Specify any requirements, such as hardware features or the operating system version, that devices must have before they can install this deployment type. Configuration Manager verifies that these requirements are met before content is deployed to the device.

Requirements:


| Requirement Type | Operator | Values                    |
|------------------|----------|---------------------------|
| Operating system | One of   | [All Windows 10 (64-bit)] |

Add...   Edit   Delete

< Previous   Next >   Summary   Cancel

12. There are no **Dependencies** to configure at this time.13. Click **Summary** to review the summary screen and click **Next**.

14. Click **Close** to close the **Completion** screen.



15. Close the remaining application creation dialog boxes.

# Technical Support

If you encounter any problems during the installation, please contact DTX Studio Technical Support:

[support@dtxstudio.com](mailto:support@dtxstudio.com)

<https://tech.dexis.com>

# Annex I: End User License Agreement (EULA)

## END USER LICENSE AGREEMENT

FOR DTX Studio™ Clinic

Version 4.10.x (2026-05-18)

THIS IS A CONTRACT. PLEASE READ THIS LICENSE AGREEMENT CAREFULLY BEFORE PRESSING THE “AGREE” BUTTON BELOW. BY PRESSING “AGREE”, YOU AGREE TO BE BOUND BY THE TERMS OF THIS LICENSE AGREEMENT AND ANY TERM INCORPORATED HEREIN. IF YOU DO NOT AGREE TO THE TERMS OF THIS LICENSE AGREEMENT, PRESS “DISAGREE” AND PROMPTLY RETURN THE SOFTWARE TO COMPANY; YOUR LICENSE FEE WILL THEN BE REFUNDED, IF ALREADY PAID.

### 1. LICENSE, RIGHT TO INSTALL AND USE OF THE SOFTWARE

This end user license agreement (the “Agreement”) is concluded between you and Medicim NV, Stationsstraat 102 b6, B-2800 Mechelen, Belgium and/or any other affiliated company of our group from which you received DTX Studio Clinic (“Company”) and governs the use of DTX Studio Clinic software products including related data and libraries that are part of the DTX Studio Clinic product, driver software for image acquisition devices that are part of the installer, and related written materials (the “Software”). This Agreement also governs the support and maintenance relating to the Software. The Company grants you as the end user a non-exclusive and limited license to use this copy of the Software (the “License”). The Software is licensed, not sold, to you. The License is non-transferable except for transfers of the Software to another computer or another operating system of the same user, practice location or Institution, subject to your specific license terms. The License is subject to further terms and conditions depending on the license type you have concluded with the Company as defined in section 2 below, and as defined and communicated to you by the Company prior to the conclusion of this Agreement.

### 2. LICENSE AND DEFINITIONS

#### 2.1 General license terms

**Use of the Software:** The copyright and other proprietary notices that appear on or in the Software must be reproduced on any permitted copy. You shall comply with all applicable laws and regulations in connection with the use of the Software. You may use the Software only for the intended purposes as defined in the Instructions for Use (“IFU”) and/or in the Software User’s Guide at each point in time, or for personal or internal training, education or demonstration. Any right to use the Software for external training, education or demonstration of the Software is subject to prior written approval by the Company. The Company may control the number and type of Licenses and the use of the Software by key codes or by other means.

**Authorized Users:** The Software is only available to authorized users, defined as dentists (including, but not limited to, periodontists, prosthodontists, orthodontists, endodontists), oral and maxillofacial surgeons, dental technicians, dental students under supervision, Universities or similar institutions, and dental laboratories (to the extent permitted by local applicable legislation), or authorized employees, auxiliary persons, contractors, representatives or agents thereof (“Authorized Users”). By installing, copying or by using the Software, you represent being an Authorized User. If you cannot adhere to the above representation, the Company is unwilling to license the Software and the installation and/or use of the Software is a violation of national and international copyright laws and conventions. If you accept the terms of this Agreement on behalf of a person or entity which is an Authorized User and on whose behalf you are authorized to act, you may only use the Software on behalf of such person or entity. In this case you are jointly and severally bound under this Agreement together with such person or entity.

**Entry Level Set:** A set of basic functional features of the Software as specified in your Purchase Offer.

**Initial License Period:** A period of time following your first purchase of a License to the Software. The duration of the Initial License Period is specified in the Purchase Offer.

**Connectivity:** The Software may require regular connection with our servers to maintain proper functionality.

For an Online License, Company may temporarily, without prior notice, suspend access to the Software or such Connectivity dependent features if Company does not detect a connection of the computer device to its servers for more than fourteen (14) consecutive days. Such suspension shall last until the computer device is connected again to Company's servers.

For an Offline License, the foregoing suspension mechanism does not apply to the core functionality licensed under the Offline License. However, any features or services that require Connectivity are not included in, and will not function under, an Offline License.

**Cloud Functionality:** The Software may comprise functionality, which enables either one of transmitting, receiving or storing data via or on a cloud platform managed by Company. For avoidance of doubt, depending on the nature of the Software, the type of license or the user account, the Cloud Functionality may not be available or be limited to one of transmitting, receiving or storing data.

**Offline License:** means a License under which the Software is licensed for use on designated computer device(s) without requiring Connectivity on such device(s) for continued core functionality.

**Online License:** means a License under which the Software's use and/or continued functionality requires Connectivity, including for purposes such as license validation, access to Cloud Functionality, delivery of Updates or Upgrades, and access to features or services that depend on Connectivity.

**Purchase Offer:** Documentation you received from Company at the time of or following your order of a License to the Software. This documentation may include but is not limited to, information on the Software, the purchase conditions, duration of the Initial License Period, specification of the features comprised in the Entry Level Set, specification of any additional Subscription Modules available to you in the Initial License Period. In case you acquired a License to the Software in combination with products of entities affiliated with Company, such as DEXIS imaging devices or Nobel Biocare products, the Purchase Offer may specify terms and conditions specific to such combined purchase.

**Service Provider:** a third-party service provider, from which certain services can be requested through the Software. Such services may include diagnostic, prosthetic design and surgery planning services, clinical tools, including custom tools, such as surgical guides, dental products, including implants and custom dental prosthetics based on data representing the patient's oral anatomy, ordered products, and CAD information.

**Software Version:** means a state of the Software labelled with a unique version number in the format of {major release}.{minor release} (e.g. 4.3). Where this Agreement mentions the two most recent Software Versions, it concerns the most current Software Version (x.y) and the immediate previous Software Version, wherein the immediate previous Software Version is the previous minor release (x.[y-1]) or when the most current version is a first major release (x.0) the last minor release of the previous major release ([x-1].z).

**Subscription Fee:** means the amount you pay for using the functionality provided in a Software Module for the duration of a specified period ("Subscription Period"). The pricing of the respective Subscription Modules and the associated Subscription Periods is available on our website or on the 'software' page of your DTX Studio Go account.

**Subscription Module:** means a specific combination of Software features and services adding onto the Entry Level Set to which you can subscribe. Every Subscription Module may have a different Subscription Fee.

**Update:** shall mean minor extensions to the supplied functionality and/or replacement of the manner of functioning of a supplied functionality by a newer functionality which has improved characteristics, but which does essentially retain the original functional scope and/or implements it differently.

**Upgrade:** shall mean extension of the Software feature set.

User Account: each Authorized User will receive a username and a password to create an account. It is the responsibility of each Authorized User to maintain the confidentiality of such credentials and to update them regularly to prevent unauthorized access. You agree to have a secure and strong password and to apply and maintain appropriate security measures to protect it. It is recommended that you have an active and up-to-date antivirus and anti-malware software with a correctly configured firewall on the computer on which you are using the Software. You are solely responsible for the access to your user account and any use of the Software through your user account. Any changes of records made by using your user account will be traced and assigned to you as the owner of such a user account.

## 2.2 Modules, features and Specific license terms

### 2.2.1 Entry Level Set

Company grants you a license on an Entry Level Set as specified in the Purchase Offer. During the Initial License Period and, for an Online License, subject to sufficient Connectivity, Company shall make available, free of charge, all Updates necessary for maintaining the features of the Entry Level Set. If upon expiration of the Initial License Period, you do not subscribe to a Subscription Module offering Updates, you will retain the right to use the Entry Level Set as available in the Software Version installed on your system at that time.

### 2.2.2 Subscription Modules

The Purchase Offer specifies the Subscription Modules available to you during the Initial License Period. After the Initial License Period, you can renew your subscription to one or more of these Subscription Modules by paying the applicable Subscription Fees. This renewal can occur upon the expiration of either the Initial License Period or any subsequent Subscription Period.

If a Subscription Module is not included in your Purchase Offer or your subscription to a Subscription Module lapsed, you can obtain or re-obtain access to such module for the duration of a Subscription Period by paying the applicable Subscription Fee.

### 2.2.3 Fee per use features

Your License may include features that are subject to the payment of a fee per use.

### 2.2.4 Specific license terms

In the Purchase Offer Company may have specified further terms and conditions of the specific license granted to you, including but not limited to number of Authorized Users and the number of computers on which the Software can be installed under the license.

### 2.2.5 Offline License (specific license terms)

(a) Grant. Where the Purchase Offer specifies that you have acquired an Offline License and subject to the terms of this Agreement, Company grants you a non exclusive and limited right to use the Software on the designated computer device(s) and for the number of Authorized Users specified in the Purchase Offer, without requiring Connectivity on such device(s) for the core functionality included in the Offline License.

(b) Exclusion of Connectivity dependent features. An Offline License does not include any features, services or functionality of the Software that require Connectivity, including, without limitation, Cloud Functionality, online ordering, synchronization, transmission, receiving or cloud based storage of data.

(c) Yearly subscription and renewal conditions. An Offline License is subject to a Subscription Fee and Subscription Period, as specified in the Purchase Offer. To maintain full functionality of the Offline License for each subsequent Subscription Period, you must (i) timely pay the applicable Subscription Fee and (ii) install the Update or Upgrade provided by Company for the relevant renewal period.

(d) Effect of non renewal or failure to install Update. Failure to timely renew your subscription by paying the applicable Subscription Fee and/or to timely install the required Update or Upgrade may result, to the extent permitted by applicable law, in limitation, suspension or loss of full functionality of the Offline License.

(e) Offline delivery of Updates. Company may make Updates or Upgrades for Offline Licenses available through delivery mechanisms that do not require Connectivity on the designated device(s), including offline installation packages or other reasonable means determined by Company.

### 2.3. Payments

#### 2.3.1 Payment by credit card

If you are paying by credit card, you authorize us to charge your credit card for all fees payable during the Subscription Period. You further authorize us to use a third party to process payments, and consent to the disclosure of your payment information to such a third party.

#### 2.3.2 Payment Information

You will keep your contact information, billing information, and credit card information (where applicable) up to date. Changes may be made on the 'Practice setup' page of your DTX Studio Go account (go.dtxstudio.com). All payment obligations are non-cancellable and all amounts paid are non-refundable, except as specifically provided for in this Agreement. All fees are due and payable at the start of a Subscription Period.

## 3. TECHNICAL SUPPORT

### 3.1 Basic Technical Support

During the Initial License Period and thereafter if you have subscribed to a Subscription Module giving you access to support services, the Company may directly or through a mandated third-party provide limited basic technical support related to the usage of the Software, which may include internet-based support during business hours, and access to certain information online. The Company collects and links data about your software usage to your license so that the company can improve the technical support it provides to you. This software usage data is used by our support teams to troubleshoot support issues that you raise. For more information about the information the Company collects and how to opt out please refer to Clause 5 Data Privacy.

Basic technical support does not include

- (i) training and educational support;
- (ii) the free providing of updates, upgrades or modified versions of the Software;
- (iii) support for Software Versions which are either (i) older than the two (2) most recent Software Versions or (ii) older than three (3) years, of which the Company does not guarantee the functionality;
- (iv) support in connection with third party software. You acknowledge that additional third-party software not approved by Company for use in connection with the Software may affect the performance of the Software.
- (v) On-site support.

### 3.2 Maintenance

During the Initial License Period and thereafter if you have subscribed to a Subscription Module giving you access to Updates, Company shall provide regular Updates and Upgrades for maintaining the Software as needed for your use of the Entry Level Set and the Subscription Modules you subscribed to. The provision of these Updates and Upgrades requires sufficient Connectivity for an Online License and, for an Offline License, requires installation of Updates or Upgrades as made available by Company for offline use, and is only available for the two most recent Software Versions or Software Versions that at the time of the provision of an Update are not older than 3 years.

#### 4. ORDERING THROUGH THE SOFTWARE

Ordering through the Software and any related services or functionality that require Connectivity are not available under an Offline License. The Software may allow you to order certain services, including diagnostic, prosthetic design and surgery planning services, clinical tools, including custom tools, such as surgical guides, dental products, including implants and custom dental prosthetics based on data representing the patient's oral anatomy, ordered products, and CAD information. Any information on products provided in the Software does not constitute an offer of or solicitation for sale or other transaction. Some products shown or described in the Software may not be available or not be regulatory cleared, released or licensed for sale in all markets. For current product assortment and availability of the Company products, please contact your local Company office or dealer. Orders for dental products or services are binding for you and become binding for the Company or Service Provider upon acceptance of the order. The Company or Service Provider may reject such orders in case the provided data or information are not in line with required specifications expectations or for other legitimate reasons.

By sending an order for custom tools or dental products you represent that: (i) you have complied with the specific technical constraints that govern the manufacturing of the specific product, (ii) you are satisfied with and approve the design of the ordered product, (iii) you acknowledge that the Company or Service Provider does not make any medical or clinical review of your design, and (iv) you agree that the Company or Service Provider therefore does not assume any liability for damages or injuries resulting from your design and produced product. In addition to the provisions of this Agreement, such orders will be governed by the terms and condition of the entity you are ordering from; whereas in case of conflicts the provisions of this Agreement shall prevail.

#### 5. DATA PRIVACY

##### 5.1 Transmission of Data through the Software

As an Authorized User, you may be subject to the laws and regulations of one or more jurisdictions, including but not limited to, laws and regulations that may apply to your use, collection, disclosure, transmission or retention of personally identifiable information, including health information. It is the Authorized User's sole responsibility to comply fully with all such applicable laws and regulations.

When sending data or information to the Company, or when allowing the Company to access such data or information for the purposes of providing technical support, you expressly agree that the Company does not desire or require access to any of your patients' personally identifiable information, including health information, and that you will not submit any patients' personally identifiable information, including health information, to the Company, or refer to any particular patient, any particular patient's medical condition, or otherwise disclose to the Company any other information which could identify a particular patient, in connection with your use of the Software or the Company's provision of technical support or maintenance of the Software. The Company reserves the right to delete without prior notice the content or parts of the content of any transmission in relation to the Software sent to the Company that the Company believes is in violation of applicable laws. The Company is entitled to log and store all data pertaining to any such transmission. To the extent any such transmission contains patient data or other personal and/or confidential data, you have the full responsibility for the processing of such data and you should comply with all applicable laws and regulations, including, but not limited to, laws and regulations that require that you obtain and maintain the prior express consent of the patient before transmitting the patient's data to third parties. In exceptional circumstances we may inform you that we are unable to resolve a support request without access to patient personally identifiable information. If it is necessary for the Company to access patient personally identifiable information in order to resolve a support request we will request your approval to do so and the information will be processed in accordance with the Data Processing Agreement and/or Business Associate Agreement (as applicable set out at Annex 1 and Annex 2). You agree to indemnify and hold the Company harmless with respect to any claim based upon the content of your transmissions or the sharing of patient records. Further, you acknowledge that the Company may use data generated through

your use of the Software, as well as data made available to the Company in the context of a support request, for the purpose understanding customer use of the Software, improving the Software, and researching and developing related products and services. Where the Company has reason to believe that the data includes personally identifiable information, it will take reasonable steps to anonymise such information prior to re-use.

The Software may enable you to store and transmit records through Company's Cloud Functionality. Records that you store or transmit using the Cloud Functionality may contain patient personally identifiable information, including health information. The information will be processed in accordance with the Data Processing Agreement and/or Business Associate Agreement (as applicable set out at Annex 1 and Annex 2). Company has put in place appropriate technical and organizational measures in place in order to ensure an proportionate level of security for your data. The servers we use are based in the European Economic Area.

If you are practicing in the United States, the servers we use are based in the United States and by accepting these terms, you agree to be bound by our business associate agreement that you can find as Annex 1 to this Agreement (HIPAA Business Associate Agreement). This business associate agreement covers our relationship to you with regards to certain protected health information, in accordance with the Health Insurance Portability and Accountability Act (HIPAA).

If you are practicing in a country of the European Union, UK, Norway, Lichtenstein, Iceland or Switzerland, the servers we use are based in the European Economic Area and by accepting these terms, you agree to be bound by our data processing agreement that you can find under Annex 2 to this Agreement (EU Data Processing Addendum). This agreement covers our relationship when Company acts as a processor of the personal data you control, in accordance with the EU General Data Protection Regulation (and for the UK, the UK GDPR & DPA 2018) or other applicable data protection law.

Patient records may contain patient data. Please note, that you alone decide with whom you want to share your patient records and it is your responsibility to determine whether you require the respective patient's consent for sharing patient data with third parties (which may include Software users of your clinic, group practice or similar institution), and to obtain all required consents. We may transfer, upon your instructions (for example if you decide to share a record with a third party), such data to another jurisdiction. It is your entire responsibility to ensure that you have the right to instruct us to transfer such data to another jurisdiction and you agree to indemnify and hold the Company harmless against any claim arising out of or related to such cross-border transfer of data that would originate from your instructions. By sharing patient records using the Software via the Cloud Functionality or otherwise, including in the context of a support request, you represent and warrant that (i) you have obtained all contained patient data lawfully and with the patient's consent, if necessary, including by providing notice to patients informing them about the processing of their personally identifiable information in the context of the Software, as well as the potential disclosure of that information to the Company for the purposes anticipated in this section (ii) you have obtained the prior approval, as required, of every patient for sharing the respective patient data, (iii) you use such patient records only for lawful purposes and within the course of dental practice, and (iv) you are compliant with applicable patient rights regulation and legislation, and data protection laws. Please note that it may be possible for you to anonymize the data prior to sharing it with others by deleting the relevant patient data or by using the anonymization function in the Software. Notwithstanding anything to the contrary, it is your sole responsibility to ensure that your deletion of data or use of the anonymization function in the Software is in compliance with laws and regulations, including requirements for anonymization or de-identification of personally identifiable health information, that are applicable to you.

## 5.2 Your privacy

The Company is committed to protecting your privacy. For the avoidance of doubt, the following section applies to our use of your personal data, and not to patient data.

Type of data collected and purpose of collection: Company may collect your e-mail address, name, address, telephone number, profession or other information in relation to purchases (e.g. credit card number, user name/ user ID, IP address, MAC address or other device identifier or other information that you provide us. Such personal data may be used to improve and operate the Software including product ordering, customer service and other services including marketing (e.g. sending product information), for software enhancement, commercial business analysis and evaluation or other similar reason and to communicate with you. In some instances, such as during the installation of the Software or as part of a customer service enquiry, you may be obliged to provide us with some information in order for us to perform our services. In other instances, where permitted by applicable laws we may use your contact details in order to send you product, or Company information that we believe is of interest to you.

If you do not want to receive such marketing communication from the Company, you may inform us by sending an email at [privacy@dtxstudio.com](mailto:privacy@dtxstudio.com)

Company may also collect certain information on the use of the Software such as technical information (e.g. CPU speed and operating system), general usability information (e.g. usage of workspaces or shortcuts), actual usage information (e.g. scanner and operating system used) and commercial information (e.g. design and ordering pattern). Company may use this information, as well as anonymised information collected in the course of providing technical support, to help us troubleshoot issues, to improve the Software and customer experience, and to research and develop related products and services. Company may also use actual usage and commercial information to help us understand the use of the Software for the purposes of improving the Software in relation to the needs of our customers and for commercial business analysis and evaluation.

You may opt out of the collection of this information by sending an email to [privacy@dtxstudio.com](mailto:privacy@dtxstudio.com).

Further information on how we handle personal data can be found in our Privacy Notice: <https://envistaco.com/en/privacy-policy>

Transfer of data to third parties: Company does not sell or disclose your Personal Data to third parties without your consent, except:

- To affiliates and third-party service providers to provide services (including support services and including the third party cloud hosting provider we use for the Cloud Functionality) in connection with this Agreement and to provide information on our Software and related products, including marketing and advertising, and to support our business operations. We require these parties to handle Personal Data in accordance with our Privacy Notice.
- To another company in connection with the sale or transfer of one of our product lines or divisions, which includes the services provided through one or more of the Company's affiliates.
- To governing regulatory authorities, including the US Food and Drug Administration, or as may otherwise be necessary for Company to comply with a legal obligation or demand.
- To authorities or third parties in order to enforce this Agreement or protect the rights, property or safety of the Company, its affiliates, employees, users of the Software or the public

In the event Company decides to reorganize or divest our business through sale, merger, or acquisition, Company may share Personal Data about you with actual or prospective purchasers. We will require any actual or prospective purchasers to treat this Personal Data in a manner consistent with this notice.

You can find all related affiliates here: <https://www.envistaco.com>

Data Transfers: Company may transfer your personal data to, or store it at, a destination outside the European Economic Area (“EEA”) Switzerland or UK. It may also be processed by staff operating outside the EEA, Switzerland or UK, who work for us or for one of our suppliers. Such staff may be engaged in, among other things, the fulfilment of your order, the processing of your payment details and the provision of support services. Where the Company has to transfer your personal data to other third countries, a legally-recognized derogation (for occasional or one-off transfers) we will use appropriate approved safeguards or we will seek your explicit consent. For further information, please do not hesitate to contact us using the details below:

Medicim NV, Stationsstraat 102, 2800 Mechelen, Belgium, or contact us at [privacy@dtxstudio.com](mailto:privacy@dtxstudio.com)

Your rights: You have the right to see and get a copy of Personal Data about you that we maintain as well as to ask us to make corrections to inaccurate or incomplete Personal Data about you. You have the right to receive data you have provided to us in a machine-readable format and to transmit that data to another controller. You may also request the erasure of your Personal Data or the restriction of its processing, object to the processing of Personal Data about you or withdraw your consent. To seek access to Personal Data about you, to file a concern, complaint, or request for correction, or to opt-out of particular programs, please contact our Privacy Office by emailing us at [privacy@dtxstudio.com](mailto:privacy@dtxstudio.com).

Your Local Data Protection Authority is responsible for making sure that privacy law is followed in the relevant countries. For more information about your privacy rights, or if you are not able to resolve a problem directly with us and wish to make a complaint, contact your local Data Protection Authority: [http://ec.europa.eu/justice/data-protection/bodies/authorities/index\\_en.htm](http://ec.europa.eu/justice/data-protection/bodies/authorities/index_en.htm).

Security Measures: The Company will implement appropriate technical and organizational measures to protect personal data against accidental or unlawful destruction or accidental loss, alteration, unauthorized disclosure or access by unauthorized third parties and against other unlawful forms of processing. However, we cannot guarantee that unauthorized third parties will not thwart our security systems and gain unlawful access to your personal data.

## 6. COPYRIGHT AND OTHER INTELLECTUAL PROPERTY

The Company and its affiliates own or have the right to license the Software and the Software is protected by national copyright law and international treaty provisions. You may not modify, adapt, translate, reverse engineer, decompile, disassemble, or otherwise attempt to discover the source code of the Software (except to the extent required to obtain interoperability with other independently created computer programs or as permitted by compulsory law). The Software may also be covered by one or more patents.

The Software comprises references to third-party products. These products are referred to using the tradenames and trademarks owned by such third parties.

## 7. RESTRICTIONS

Except as expressly permitted by this Agreement or by the applicable mandatory law, you are not entitled to lease, rent, loan, timeshare, resell, assign, transfer, license or sublicense, or otherwise distribute or grant any rights of use (such as usufruct) to the License or the Software to any third party without the explicit written consent of the Company.

You agree that the Software will not be shipped, transferred or exported into any country or used in any manner prohibited by the United States Export Administration Act, or any other export laws, restrictions or regulations.

These items are controlled and authorized by the U.S. Government for export only to the specified country of ultimate destination for use by the end-user herein identified. They may not be resold, transferred, or otherwise disposed of, to any other country or to any person other than the authorized end-user or consignee(s), either in their original form or after being incorporated into other items, without first obtaining approval from the U.S. government or as otherwise authorized by U.S. law and regulations.

## 8. WARRANTY AND LIABILITY

The Company warrants, for a period of six (6) months from delivery, that the Software as made available to you by Company will perform in accordance with the Company's then current product description when operated according to its instructions. The Company's liability for any damages caused by the use of the Software is excluded to the extent permitted by law, if you have omitted to Update/Upgrade your installed copy of the Software for six (6) months or more with such Update/Upgrade as the Company has made available in accordance with section 3 above. Further Company's liability under this Agreement is limited to an obligation of the Company to correct, or in its sole discretion, to replace the Software with software of equivalent functionality. Correction and replacement, as specified in this section 8, are the sole remedies available for you in case of any defects in the Software.

Company will and cannot guarantee the functionality of a Software Versions (i) older than the two (2) most recent Software Versions or (ii) older than three (3) years.

By installing and/or using the software, you confirm that you have the full knowledge how to use the Software safely and effectively. The Software is no replacement of the adequate and careful dental diagnostics and treatment in accordance with the medical duty of care. You are using the Software at your own risk. The company therefore specifically denies any liability for indirect, incidental, special or consequential damages resulting from the sale or use of the software, except if such damage is caused by a gross negligent or wilful misconduct of the company to the extent permitted by applicable law.

## 9. INDEMNIFICATION

The Company will defend, indemnify, save and hold you harmless from any money judgment, costs, and attorney fees awarded or in settlement to the extent such are due to a claim that the Software, as provided by the Company, infringes a copyright or protected trade secret of any third party. This indemnity is subject to the Company being provided prompt notice of such claim and being granted full control of the defence or settlement thereof. The Company specifically disclaims any liability for patent infringement. The indemnity of this section 9 is the sole remedy available for you for infringement of any rights.

You shall indemnify and hold harmless the Company from any damages and any money judgment, costs, and attorney fees awarded or in settlement to the extent such are due to a claim from a third party based on the unlawful disclosure, use, sharing or transmission of patient data or the infringement or breach of applicable patient rights regulation and legislation, and data protection laws, or for any breach of this Agreement by you.

## 10. TERM AND TERMINATION

This Agreement and the herewith agreed upon License are effective from the date of your acceptance (which is submitted by pressing the "AGREE" button below) until its termination in accordance with this section 10.

With respect to an Offline License, failure to timely pay the yearly subscription fee and/or failure to install the Update or Upgrade required for the applicable renewal period, as set out in Section 2.2.5, may result in suspension or limitation of full functionality and may constitute non-compliance with this Agreement for purposes of this Section 10.

Right to terminate: Company or you may terminate the Agreement and the License unilaterally at any time by prior written notice to the other Party, including e-mail or other electronic communication in connection with the Software.

Automatic termination: The License shall automatically terminate upon failure by you to comply with any term of this Agreement.

Upon termination: you agree to immediately destroy the Software and all copies, modifications and merged portions in any form and promptly return the Software to the Company. Your User Account will be deleted promptly. The Company shall block any Authorized User from accessing the Software under the same License by any appropriate technical means. The Company will store records received from you for ninety (90) days from the termination in which period you shall be entitled to obtain a copy of any such data or records. After this period, the Company will delete such records, unless required otherwise.

#### 11. GENERAL PROVISIONS, GOVERNING LAW AND JURISDICTION

If any part of this Agreement is found void and unenforceable, the remaining of the Agreement shall remain valid and enforceable according to its terms. The Company is entitled to change the terms of this Agreement; any and all changes shall become effective the earliest of You accepting them (e.g. at installation of an Update/Upgrade), or thirty (30) days after the communication of such change (by notice to you or by update of these terms). This Agreement shall be governed by the substantive laws of Belgium. The exclusive place of jurisdiction shall be at the registered office of Company.

#### Annex 1: HIPAA Business Associate Agreement

This Business Associate Agreement (“BAA”) is entered into by and between the authorized users (as defined in the End User License Agreement) (“Covered Entity”) and Medicim NV, Stationsstraat 102 b6, B-2800 Mechelen, Belgium and/or any other affiliated company of our group from which you received DTX Studio Clinic (collectively “Business Associate”) (collectively referred to herein as the “Parties”), effective as of the date the End User License Agreement is agreed upon by the Covered Entity.

WHEREAS, the Parties have entered into one or more services agreements (“Services Agreement”) whereby Business Associate performs certain functions, activities, or services for or on behalf of Covered Entity that involve the use or disclosure of Protected Health Information (as defined herein) and Electronic Protected Health Information (as defined herein); and

WHEREAS, this BAA is intended to comply with the requirements of the Health Insurance Portability and Accountability Act of 1996 (“HIPAA”) and implementing regulations, the Standards for Privacy of Individually Identifiable Health Information (the “Privacy Rule”) the Security Standards for the Protection of Electronic Protected Health Information (the “Security Rule”), and the privacy, security and Breach Notification regulations of the Health Information Technology for Economic and Clinical Health (“HITECH”) Act and the HIPAA Omnibus final rule (collectively, the “HIPAA Rules”), as amended from time to time.

NOW, THEREFORE, in consideration of the Parties’ continuing obligations under the Services Agreement between the Parties, the agreements herein, and other good and valuable consideration, the receipt and sufficiency of which is hereby acknowledged, the Parties agree as follows:

##### 1. Definitions

Except as otherwise defined herein, any and all capitalized terms in this BAA shall have the definitions set forth in the Privacy Rule or the Security Rule.

- (A) “Breach” has the meaning given to such term in 45 C.F.R. § 164.402.
- (B) “Business Associate” has the meaning set forth above.
- (C) “Covered Entity” has the meaning set forth above.
- (D) “Designated Record Set” has the same meaning as the term “designated record set” in 45 C.F.R. § 164.501 of the Privacy Rule.
- (E) “Electronic Protected Health Information” (“EPHI”) has the same meaning as the term “electronic protected health information” in 45 C.F.R. § 160.103 of the Security Rule, limited to the information created or received by Business Associate from or on behalf of Covered Entity.

- (F) “Health Information Technology for Economic and Clinical Health (“HITECH”) Act” has the meaning set forth above.
- (G) “HIPAA” has the meaning set forth above
- (H) “Individual” has the same meaning as the term “individual” in 45 C.F.R. § 160.103 of the Privacy Rule.
- (I) “Privacy Rule” has the meaning set forth above.
- (J) “Protected Health Information (“PHI”)” has the same meaning as the term “protected health information” in 45 C.F.R. § 160.103 of the Privacy Rule (including, without limitation, Electronic Protected Health Information), limited to the information created or received by Business Associate from or on behalf of Covered Entity.
- (K) “Required by Law” has the same meaning as the term “required by law” in 45 C.F.R. § 164.103 of the Privacy Rule.
- (L) “Secretary” means the Secretary of the Department of Health and Human Services or his or her designee.
- (M) “Security Incident” has the same meaning as the term “security incident” in 45 C.F.R. § 164.304 of the Security Rule.
- (N) “Security Rule” has the meaning set forth above.
- (O) “Unsecured PHI” has the meaning given to such phrase in the Breach Notification Rule at 45 C.F.R. § 164.402.

## 2. Obligations and Activities of Business Associate

- (A) Business Associate acknowledges and agrees that all PHI that is created or received by Covered Entity and used by or disclosed to Business Associate or created or received by Business Associate on Covered Entity’s behalf shall be subject to this BAA.
- (B) Business Associate agrees to not use or disclose PHI other than as permitted or required by this BAA or as Required by Law.
- (C) Business Associate agrees to use appropriate safeguards to prevent use or disclosure of PHI other than as provided for by this BAA.
- (D) Business Associate agrees to notify Covered Entity promptly following discovery of any Breach of Unsecured PHI. Business Associate will provide such information to Covered Entity as required in the Breach Notification Rule.
- (E) Business Associate agrees to report to Covered Entity any use or disclosure of PHI not provided for by this BAA or any Security Incident of which it becomes aware.
- (F) Business Associate agrees to ensure that any agent, including a subcontractor, to whom it provides PHI received from, or created or received by Business Associate for, or on behalf of, Covered Entity agrees in writing to substantially similar restrictions and conditions that apply through this BAA to Business Associate with respect to such information.
- (G) To the extent Business Associate maintains PHI in a Designated Record Set, Business Associate will make such PHI available to Covered Entity within thirty (30) business days of a request by Covered Entity for access to such PHI. For avoidance of doubt, Covered Entity understands and agrees that Business Associate does not maintain any PHI in a Designated Record Set. If an Individual makes a request for access directly to Business Associate, Business Associate will within thirty (30) business days forward such request in writing to Covered Entity. Covered Entity will be responsible for making all determinations regarding the grant or denial of an Individual’s request for PHI and Business Associate will make no such determinations. Only Covered Entity will release PHI to an Individual pursuant to such a request, unless Covered Entity directs Business Associate to do so.

(H) To the extent Business Associate maintains PHI in a Designated Record Set, Business Associate will provide such PHI to Covered Entity for amendment within thirty (30) business days of receiving a request from Covered Entity to amend an Individual's PHI. For avoidance of doubt, Covered Entity understands and agrees that Business Associate does not maintain any PHI in a Designated Record Set. If an Individual makes a request for amendment directly to Business Associate, Business Associate will within thirty (30) business days forward such request in writing to Covered Entity. Covered Entity will be responsible for making all determinations regarding amendments to PHI and Business Associate will make no such determinations unless Covered Entity directs Business Associate to do so.

(I) Within thirty (30) days of receiving a written request from Covered Entity, Business Associate shall provide to Covered Entity an accounting of the disclosures of the Individual's PHI in accordance with 45 C.F.R. § 164.528. If an Individual requests an accounting of Disclosures directly from Business Associate, Business Associate will forward the request and its record of Disclosures to Covered Entity within thirty (30) business days of Business Associate's receipt of the Individual's request. Covered Entity will be responsible for preparing and delivering the accounting to the Individual. Business Associate will not provide an accounting of its Disclosures directly to any Individual, unless directed by Covered Entity to do so.

(J) Business Associate agrees to make internal practices, books, and records, including policies and procedures and PHI, relating to the use and disclosure of PHI received from, or created or received by Business Associate on behalf of, Covered Entity, available to the Secretary, in a time and manner designated by the Secretary, for purposes of the Secretary determining Covered Entity's compliance with the Privacy Rule.

### 3. Permitted Uses and Disclosures by Business Associate

(A) Except as otherwise limited by this BAA, Business Associate may use or disclose PHI to perform functions, activities or services for or on behalf of Covered Entity as contemplated in the Services BAA, provided that such use or disclosure does not violate the Privacy Rule or the HITECH Act if done by Covered Entity.

(B) Except as otherwise limited by this BAA, Business Associate may use PHI for the proper management and administration of Business Associate or to carry out the present and/or future legal responsibilities of Business Associate.

(C) Except as otherwise limited by this BAA, Business Associate may disclose PHI for the proper management and administration of Business Associate, provided that disclosures are Required by Law, or Business Associate obtains reasonable assurances from the person to whom the PHI is disclosed that it will remain confidential and be used or further disclosed only as Required by Law or for the purpose for which it was disclosed to the person, and the person notifies Business Associate of any breaches in the confidentiality of the PHI.

(D) Business Associate may use PHI to report violations of law or other conduct to appropriate federal and state authorities or other designated officials, consistent with 45 C.F.R. § 164.502(j)(1).

(E) Business Associate may use PHI to aggregate data as permitted by 45 C.F.R. § 164.504(e)(2)(i)(B).

(F) Business Associate may use PHI to create de-identified information in accordance with 45 CFR § 164.514.

### 4. Obligations of Covered Entity on Behalf of Business Associate

(A) Covered Entity shall notify Business Associate of any limitation(s) in its notice of privacy practices within fifteen (15) business days of Covered Entity's receipt of the Individual's request in accordance with 45 C.F.R. § 164.520, to the extent that such limitation(s) may affect Business Associate's use or disclosure of PHI.

(B) Covered Entity shall notify Business Associate of any changes in, or revocation of, permission by an Individual to use or disclose PHI within fifteen (15) business days of Covered Entity's receipt of the Individual's request, to the extent that such changes may affect Business Associate's use or disclosure of PHI.

(C) Covered Entity shall notify Business Associate of any restriction to the use or disclosure of PHI that it has agreed to within fifteen (15) business days of Covered Entity agreeing to such restriction in accordance with 45 C.F.R. § 164.522, to the extent that such restriction may affect Business Associate's use or disclosure of PHI.

(D) Covered Entity shall not request that Business Associate use or disclose PHI in any manner that would not be permissible under HIPAA if done by a Covered Entity (unless permitted by HIPAA for a business associate).

(E) Covered Entity is responsible for implementing appropriate privacy and security safeguards to protect its PHI in compliance with HIPAA.

5. Security Rule and HITECH Act Responsibilities of the Business Associate.

With regard to its use and/or disclosure of ePHI, Business Associate hereby agrees to do the following:

(A) Comply with the applicable requirements of the Security Rule.

(B) Require all of its subcontractors and agents that create, receive, maintain, or transmit ePHI on behalf of Business Associate to agree, in writing, to adhere to substantially similar restrictions and conditions concerning ePHI that apply to Business Associate pursuant to Section 5 of this BAA.

(C) Report to Covered Entity any Security Incident of which it becomes aware. Specifically, Business Associate will report to Covered Entity any successful unauthorized access, Use, Disclosure, modification, or destruction of ePHI or interference with system operations in an information system containing ePHI of which Business Associate becomes aware within thirty (30) business days of Business Associate learning of such Security Incident. The parties agree that this Section serves as notice by Business Associate to Covered Entity of the ongoing existence of attempted but Unsuccessful Security Incidents (as defined below), for which no additional reporting is required. For purposes of this BAA, "Unsuccessful Security Incidents" include but are not limited to activity such as "pings" and other broadcast attacks on Business Associate's firewall, port scans, unsuccessful log-on attempts, denials of service and any other attempts to penetrate such computer networks or systems that do not result in unauthorized access, use or disclosure of ePHI.

6. Term and Termination

(A) Term. The Term of this BAA shall in effect as of the Effective Date set forth above, and shall terminate when all the PHI provided by Covered Entity to Business Associate, or created or received by Business Associate for or on behalf of Covered Entity, is destroyed or returned to Covered Entity or, if it is infeasible to return or destroy the PHI, protections are extended to such information, in accordance with the termination provisions in this Section 6.

(B) Termination for Cause. If Covered Entity or Business Associate knows of a material breach or violation by the other party of any provision of this BAA, then the non-breaching party shall provide written notice of the breach or violation to the other party that specifies the nature of the breach or violation. The breaching party must cure the breach or end the violation within thirty (30) days after receipt of the written notice. In the absence of a cure reasonably satisfactory to the non-breaching party, then the non-breaching party may terminate this BAA between the parties.

(C) Effect of Termination.

(i) Except as provided in paragraph (ii) of this Section, upon termination of this BAA, for any reason, Business Associate shall return or destroy all PHI received from Covered Entity, or created or received by Business Associate for or on behalf of Covered Entity. This provision shall apply to PHI that is in the possession of subcontractors or agents of Business Associate. Business Associate shall retain no copies of the PHI.

(ii) In the event that Business Associate determines that returning or destroying the PHI is infeasible, Business Associate shall extend the protections of this BAA to such PHI and limit further uses and disclosures of such PHI to those purposes that make the return or destruction infeasible, for so long as Business Associate maintains such PHI.

#### 7. Notification

With respect to notices pursuant to paragraph 2(D) above, notice shall be made by telephone to the telephone number associated with Covered Entity's account, followed promptly by a written notice as described below.

Any notices required or provided for under this BAA shall be made in writing and shall be either personally delivered, mailed by first class mail or sent via facsimile or electronic mail to the appropriate individual identified below:

For Covered Entity:

Your address

For Business Associate:

Medicim NV, Stationsstraat 102 b6, B-2800 Mechelen, Belgium or contact us at [privacy@dtxstudio.com](mailto:privacy@dtxstudio.com)

Either Party may designate a different address in writing to the other.

#### 8. Regulatory References

A reference in this BAA to a section in the Privacy Rule, the Security Rule or the HITECH Act means the section as in effect or as amended.

#### 9. Survival

The respective rights and obligations of the Business Associate under Section 6 of this BAA shall survive the termination of this BAA.

#### 10. Interpretation

Any ambiguity in this BAA shall be resolved to permit compliance with the HIPAA Rules. Any conflict between the terms of this BAA and any other BAA relating to the same subject matter shall be resolved so that the terms of this BAA supersede and replace the relevant terms of any such other BAA.

#### 11. Counterparts

This BAA may be executed in counterparts which, when all signatures are assembled, shall have the same effect as a single, fully-executed agreement. Facsimile and photocopy signatures shall have the same binding effect as manual signatures.

#### 12. Severability

The provisions of this BAA shall be severable, and if any provision of this BAA shall be held or declared to be illegal, invalid or unenforceable, the remainder of this BAA shall continue in full force and effect as though such illegal, invalid or unenforceable provision had not been contained herein.

#### 13. Governing Law

Except to the extent that the HIPAA Rules or other federal law applies, this BAA and the obligations of the Parties hereunder will be governed by interpreted in accordance with the laws of the State of California.

14. Effect

This Agreement amends, restates and replaces in its entirety any prior business associate agreement between the parties. This Agreement supersedes all prior or contemporaneous written or oral contracts or understandings between the parties relating to their compliance with health information confidentiality laws and regulations, including HIPAA and HITECH.

15. No Agency Relationship.

It is not intended that an agency relationship (as defined under the federal common law of agency) be established hereby expressly or by implication between Covered Entity and Business Associate under HIPAA or the Privacy Rule, Security Rule, or Breach Notification Rule. No terms or conditions contained in this BAA shall be construed to make or render Business Associate an agent of Covered Entity.

ANNEX 2: EU Data Processing Addendum

This EU Data Processing Addendum ("Addendum") is made by and between Authorized User(s) to whom the Software is Licensed ("Customer") and Medicim NV, Stationsstraat 102, 2800 Mechelen, Belgium, and/or any other affiliated company of our group from which you received the Software ("Company"), and amends and is incorporated into the End User License Agreement between Customer and Company ("Agreement"). For good and valuable consideration, the receipt and sufficiency of which are hereby acknowledged, Customer and Company agree as follows:

1. This Addendum shall apply to all Processing of Personal Data performed by Company on behalf of Customer under the Agreement. In case of any direct conflict between this Addendum and the Agreement, the terms of this Addendum shall prevail. For the avoidance of doubt, this Addendum does not apply to Company's Processing of Authorized User Personal Data under Section 5.2 of the Agreement, which Company Processes for its own purposes as a Controller.
2. For the purpose of this Addendum, Controller, Processor, Data Subject, Personal Data, Personal Data Breach, and Processing have the meanings ascribed to them by Applicable Data Protection Law. Applicable Data Protection Law means all applicable UK, Swiss, EU or EU Member State laws and regulations relating to the privacy, confidentiality, security or protection of Personal Data, including, without limitation, (i) Regulation (EU) 2016/679 ("GDPR") and EU Member State laws supplementing the GDPR, (ii) the EU Directive 2002/58/EC (e-Privacy Directive), as replaced from time to time, and EU Member State laws implementing the e-Privacy Directive, including laws regulating the use of cookies and other tracking technologies, (iii) the UK GDPR, as that term is defined by section 3(10) of the Data Protection Act 2018, and (iv) the Swiss Federal Act on Data Protection of 2020.
3. In circumstances in which Company Processes Personal Data as a Data Processor under the Agreement, Company shall:
  - (a) Process the Personal Data only in accordance with the documented instructions of Customer, unless Company is required to do otherwise by applicable law and permitted to do otherwise by Applicable Data Protection Law, in which case Company shall inform Customer of the relevant legal requirement before Processing the Personal Data unless informing Customer is prohibited by applicable law on important grounds of public interest;
  - (b) Ensure that Company's personnel authorized to Process the Personal Data have committed themselves to confidentiality or are under an appropriate statutory obligation of confidentiality;
  - (c) Implement appropriate technical and organizational measures to ensure a level of security appropriate to the risk of the Processing of the Personal Data and in line with the requirements of Applicable Data Protection Law;

- (d) Taking into account the nature of the Processing, assist Customer, by appropriate technical and organizational measures, insofar as this is possible, in fulfilling Customer's obligation to respond to Data Subjects' requests for exercising their rights under the GDPR with respect to their Personal Data;
- (e) Assist Customer in complying with its obligations relating to data security, personal data breaches and data protection impact assessments, taking into account the nature of the Processing and the information available to Company;
- (f) At Customer's choice, delete or return all Personal Data to Customer after the end of the term of the Agreement, and delete existing copies, unless Applicable Law requires storage of the Personal Data;
- (g) Make available to Customer, for inspection on Company's premises only, the information necessary to demonstrate compliance with the obligations set out in this Addendum and allow for and contribute to audits conducted by Customer or another auditor mandated by Customer and approved by Company, provided that Customer gives Company at least 30 days' prior written notice of its intention to carry out an audit. This notice shall include a detailed work plan for the audit. Any third party involved in the audit must agree to Company's confidentiality undertakings and Customer will bear all costs and expenses incurred by Company in connection with the audit; and
- (h) Company shall immediately inform Customer if, in Company's opinion, an instruction provided by Customer infringes Data Protection Applicable Law.
4. Customer agrees that Company may subcontract its Processing operations performed on behalf of Customer under the Agreement. Prior to permitting any subcontractor to Process Personal Data ("Sub-Processor"), Company shall require each such Sub-Processor to enter into a written agreement that imposes the same data protection obligations as set out in this Addendum and shall remain liable to Customer for the performance of the Sub-Processor (subject to the limitations on liability set out in this Agreement). As at the date of this Agreement, the Sub-Processors used by Company are set out in Exhibit A (Description of Processing) to this Addendum. In the event that Company wishes to appoint a new Sub-Processor, it shall notify Customer (via the Software portal and/or by email) and provide Customer with an opportunity to object to the appointment.
5. For the purposes of its Processing instructions under paragraph 3(a), the Customer expressly authorises the Company to anonymise Personal Data so that the Company can process the anonymised data for its own purposes, including without limitation for improving the Software and researching and developing related products and services. The Company shall be responsible for ensuring the lawfulness of its anonymisation and its subsequent use of anonymised data.
6. Customer agrees that Company may transfer Personal Data outside of the EEA, the United Kingdom or Switzerland for the purpose of fulfilling its obligations to Customer under the Agreement and on the condition that Company has implemented appropriate safeguards for the transfer of the Personal Data in accordance with applicable law.

## Exhibit A (Description of Processing)

|                                                  |                                                                                                                                                                             |
|--------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Subject matter, nature and purpose of Processing | The transmitting, receiving or storing of data as part of the Cloud Functionality                                                                                           |
| Categories of Personal Data                      | Dental scans and other medical images; medical notes; name; patient ID number                                                                                               |
| Categories of Data Subject                       | Patients of the Customer                                                                                                                                                    |
| Duration of the Processing                       | For the duration of the Agreement                                                                                                                                           |
| Sub-Processors                                   | Microsoft Azure (cloud hosting for the Cloud Storage Services and for data sharing by the Customer with third parties, such as patients and/or other medical professionals) |



Nobel Biocare AB  
Box 5190, 402 26  
Västra Hamngatan 1,  
411 17 Göteborg,  
Sweden

[www.nobelbiocare.com](http://www.nobelbiocare.com)

**Distributed in the United States by:**

Nobel Biocare USA, LLC  
22715 Savi Ranch Parkway  
Yorba Linda, CA, 92887 USA

Rx Only



For prescription use only

Medical device



Caution



[ifu.dtxstudio.com/symbolglossary](http://ifu.dtxstudio.com/symbolglossary)  
[ifu.dtxstudio.com](http://ifu.dtxstudio.com)

Link to Online Symbol  
Glossary and IFU Portal

**For Prescription Use Only. Caution** Federal (United States) law restricts this device to sale by or on the order of a clinician, medical professional or physician.