



DTX Studio™ Implant

Version 3.6

Installation Guide

Table of Contents

Introduction	3
Disclaimer of Liability	3
Cybersecurity	3
Preparing Your Computer	3
 Where to Find Installer and Install Key	 4
 Install DTX Studio Implant on a Windows® System	 5
Install the Practice-Setup License	5
Configure the Software	5
 Install DTX Studio Implant on a Mac OS® X System	 7
Install the Practice-Setup License	7
Configure the Software	7
 Technical Support	 9

Introduction

Disclaimer of Liability

This product is part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Non-recommended use of products made by third parties in conjunction with Nobel Biocare products will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare products has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The user is also obliged to study the latest developments in regard to this Nobel Biocare product and its applications regularly. In cases of doubt, the user has to contact Nobel Biocare. Since the utilization of this product is under the control of the user, it is his/her responsibility. Nobel Biocare does not assume any liability whatsoever for damage arising thereof. Please note that some products detailed in this Instructions for Use may not be regulatory cleared, released or licensed for sale in all markets.

Cybersecurity

It is recommended that active and up-to-date antivirus and anti-malware software, together with a correctly configured firewall, are installed on the computer where DTX Studio Implant is to be used.

Furthermore, always lock your computer when it is left unattended. Failure to do so may lead to unintended handling of the diagnosis and planning or treatment.

Preparing Your Computer

Make sure your computer is connected to the Internet.

Check that your computer can open and display .pdf files (double-click on a .pdf file to test). Adobe® Acrobat® Reader is an example of software typically used to open and display .pdf files.

Log on to your computer with a user account that has administrator rights. We recommend that all other software on your computer is closed before DTX Studio Clinic is installed.

Ensure your computer meets the recommended system requirements. For detailed information refer to the DTX Studio Implant Computer Guidelines. Download these guidelines from the [DTX Studio Documentation Library](#) or [DTX Studio Go](#). If you need a digital or printed copy of the latest version, please contact your local representative.


Where to Find Installer and Install Key

1. Go to go.dtxstudio.com . The portal is best viewed with Chrome, Safari and Edge browsers. Internet Explorer is currently not supported.
2. Log in with the username and password that were sent to you. If this is the first time you log in to DTX Studio Go, you will be prompted to change the password.


Note

If you did not receive username and password by mail, please contact customer support.

3. On the [Software](#) tab, click the DTX Studio Implant version you need.
4. To download the latest version of the selected software product, click [Mac](#) or [Windows](#), depending on the operating system you wish to install DTX Studio Implant on. The installer file is downloaded automatically to your download folder.



5. Copy the installation key.



6. In this guide, go to the section of the software version you wish to install and follow the instructions:
 - “Install DTX Studio Implant on a Windows® System” on page 5
 - “Install DTX Studio Implant on a Mac OS® X System” on page 7

Install DTX Studio Implant on a Windows® System

Please follow the instructions carefully.

Have the DTX Studio Implant license information document at hand.

Note

For this license, a maximum of ten users can be created. This license can be installed on six devices.

Install the Practice-Setup License

1. Double-click DTXImplant.exe.

Note

Before starting the actual installation, the system will check whether a number of prerequisite programs are available on the system you are using. If not available yet, a message is shown and the installation of the program(s) is requested. Accept to have them installed before starting the actual installation of DTX Studio Implant.

2. A welcome message for the install wizard appears. Click **Next**.
3. The End-User License Agreement is displayed. Read it carefully. You can also print this information. To accept the agreement, click **I accept the terms in the license agreement**, and click **Next**. Alternatively, click **Cancel** and the installation procedure will abort.
4. Choose a destination folder where the DTX Studio Implant files will be stored. By default, this folder will be C:\Program Files\DTX Studio Implant. To change this destination folder, click **Change**, browse to the folder of your choice and click **OK**.
5. Click **Next**.
6. Click **Install** to start the installation.
7. When the software is installed, click **Finish**. DTX Studio Implant starts and a configuration wizard is shown, allowing you to set up the software.


Configure the Software

1. In the configuration wizard, click **Next**.
2. Enter the installation key that you copied from DTX Studio Go.
3. Click **Next**.
4. Select the team folder and the working folder. In the **Team folder** all patient records will be stored. It should be accessible to all practice members. You can still change this setting by selecting **General Preferences** from the **DTX Studio Implant** menu. In the **Working folder** the privately stored patient-related data will be saved. The default location is C:\DTX Studio\Implant. The location for this working folder cannot be changed in the **General Preferences**. Please make sure that this folder is stored on a hard disk with sufficient free disk space (approximately 250 MB per patient file is needed). All users of the computer should have read/write access to this folder.
5. Click **Next**.

Note

If there are Mac® and Windows® users in your practice, the team folder has to be set once on a Mac® computer for all Mac® users and once on a Windows® computer for all Windows® users.

6. Decide whether you would like to install products from a local file or download them from the server (a stable, high-speed Internet connection is required).
 - If you choose the local file, browse to the implant library file (with the extension .ilf) and click **Open** to install all products in the library. The installation progress is shown.
 - If you choose to download the products from the server, select **Install through the download server** and click **Next**. A new wizard page allows you to select the manufacturer(s) whose products are to be downloaded. Select the manufacturer(s) and click **Next** to start downloading. The download progress is shown. It is possible to cancel the process by returning to the previous page and switching to another option.
7. When all products are installed, click **Finish**.
8. Log on as a specific user.
9. The End-User License Agreement is displayed. Read it carefully. You can also print this information. To accept the agreement, click **I accept the terms in the license agreement**, and click **Next**. Alternatively, click **Cancel** and the installation procedure will abort.
10. When the account is successfully registered, click **Finish**. The Instructions for Use website is displayed. The software is ready to use.

A DTX Studio Implant icon  is available on your desktop. Double-click it to start the software. DTX Studio Implant can also be started by clicking the Windows® button. Browse the list of all programs that appears and click **Nobel Biocare**. Click **DTX Studio Implant** to start the software. The Instructions for Use are also available here in a printable format.

At all times, the License Agreement is available in the **Manage account** dialog within DTX Studio Implant. It is advised to use the software while connected to the Internet in order to receive software updates etc. An Internet connection is required to prepare and send orders.

Install DTX Studio Implant on a Mac OS® X System

Please follow the instructions carefully.

Have the DTX Studio Implant license information document at hand.

Note

For this license, a maximum of ten users can be created. This license can be installed on six devices.

Install the Practice-Setup License

1. Double-click the Apple® Disk Image file (.dmg) to mount it. The mounted **DTX Studio Implant** disk will appear on your desktop.
2. Double-click this file. An installation dialog appears, showing different options:
 - **Install DTX Studio Implant**
 - **Software License Agreement**
3. To install DTX Studio Implant, double-click **Install DTX Studio Implant**. A welcome screen is displayed. Click **Continue**.
4. The End-User License Agreement is displayed. Read it carefully. You can also print this information.
5. Click **Continue**. A warning is displayed, asking for acceptance of the Software License Agreement.
6. Click **Agree** to continue with the installation. If you click **Disagree** the installation procedure is aborted.
7. Keep the default location for the installation (**Applications**) or change it to the location of your choice. This is the disk where the software files will be saved. Click **Continue**.
8. Click **Install**.
9. Enter the username and password of an administrator account.
10. When the software is installed, click **Finish**. DTX Studio Implant starts and a configuration wizard is shown, allowing you to set up the software.

Configure the Software

1. In the configuration wizard, click **Continue**.
2. Enter the installation key that is printed on the license information sheet. Click **Continue**.
3. Configure the software. Select the team folder and the working folder. In the **Team folder**, all patient records will be stored. It should be accessible to all practice members. You can still change this setting by selecting **General Preferences** from the **DTX Studio Implant** menu. In the **Working folder**, the privately stored patient-related data will be saved. The default location is \\Users\\Shared\\DTX Studio Implant. The location for this folder cannot be changed in the **General Preferences**. Please make sure that this folder is stored on a hard disk with enough free disk space (approximately 250 MB per patient file is needed). All users of the computer should have read/write access to this folder.

Note

If there are Mac® and Windows® users in your practice, the team folder has to be set once on a Mac® computer for all Mac® users and once on a Windows® computer for all Windows® users.

4. Click **Continue**.

5. Decide whether you would like to install products from a local file or download them from the server (a stable, high-speed Internet connection is required).
 - If you choose the local file, browse to the implant library file (with the extension .ilf) and click **Open** to install all products in the library. The installation progress is shown.
 - Should you choose to download the products from the server, a new wizard page allows you to select the manufacturer(s) whose products are to be downloaded. Click **Continue** to start downloading. The download progress is shown. It is possible to cancel the process by returning to the previous page and switching to another option.
6. When all products are installed, click **Finish**.
7. Log on as a specific user.
8. The End-User License Agreement is displayed. Read it carefully. You can also print this information. To accept the agreement, click **I accept the terms in the license agreement** and click **Next**. Alternatively, click **Exit** and the installation procedure will abort.
9. Click **Done**. The Instructions for Use website is displayed. The software is ready to use.

A DTX Studio Implant icon  is available in **Applications** (default location). Click it to start the software.

At all times, the License Agreement is available in the **Manage account** dialog within DTX Studio Implant. It is advised to use the software while connected to the Internet in order to receive software updates etc. An Internet connection is required to prepare and send orders.

Technical Support

If you encounter any problems during the installation, please contact DTX Studio Technical Support:

- support@dtxstudio.com
- <https://tech.dexis.com/support>

For North America:

- dtxsupport.na@dexis.com
- Phone number: +1-833-389-2255



Nobel Biocare AB
Box 5190, 402 26
Västra Hamngatan 1,
411 17 Göteborg,
Sweden

www.nobelbiocare.com

Sponsor and Distributed in Australia by:

Nobel Biocare Australia Pty Ltd
Suite 4.02, Level 4, Building A,
1 Eden park drive Macquarie Park NSW 2113
Australia

Phone: +61 1800 804 597

Distributed in New Zealand by:


Nobel Biocare New Zealand Ltd
33 Spartan Road
Takanini, Auckland, 2105
New Zealand

Phone: +64 0800 441 657

Distributed in Turkey by:


EOT Dental
Sağlık Ürünleri ve Dış Ticaret A.Ş.
Nispetiye Mah. Aytar Cad.
Metro İş Merkezi No: 10/7
Beşiktaş İSTANBUL

Phone: +90 2123614901



CH importer/representative:

Nobel Biocare Services AG
Balz Zimmermann-Strasse 7
8302 Kloten
Switzerland



UK responsible person:

Nobel Biocare UK Ltd.
4 Longwalk Road
Stockley Park
Uxbridge UB11 1FE
United Kingdom



ifu.dtxstudio.com/symbolglossary
ifu.dtxstudio.com